



PRTFDC1 blocking peptide (CDBP5961)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PRTFDC1 phosphoribosyl transferase domain containing 1 [Homo sapiens (human)]
Official Symbol	PRTFDC1
Synonyms	PRTFDC1; phosphoribosyl transferase domain containing 1; HHGP; phosphoribosyltransferase domain-containing protein 1
Entrez Gene ID	56952
mRNA Refseq	NM_001282786
Protein Refseq	NP_001269715
UniProt ID	Q9NRG1
Function	NOT hypoxanthine phosphoribosyltransferase activity; NOT magnesium ion binding;

