



PTCHD2 blocking peptide (CDBP5963)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PTCHD2 patched domain containing 2 [Homo sapiens (human)]
Official Symbol	PTCHD2
Synonyms	PTCHD2; patched domain containing 2; DISP3; patched domain-containing protein 2; dispatched 3; protein dispatched homolog 3
Entrez Gene ID	57540
mRNA Refseq	NM_020780
Protein Refseq	NP_065831
UniProt ID	Q9P2K9
Function	hedgehog receptor activity; molecular_function