



# PTCHD2 blocking peptide (CDBP5964)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

## GENE INFORMATION

|                 |                                                                                                                             |
|-----------------|-----------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">PTCHD2 patched domain containing 2 [ Homo sapiens (human) ]</a>                                                 |
| Official Symbol | PTCHD2                                                                                                                      |
| Synonyms        | PTCHD2; patched domain containing 2; DISP3; patched domain-containing protein 2; dispatched 3; protein dispatched homolog 3 |
| Entrez Gene ID  | <a href="#">57540</a>                                                                                                       |
| mRNA Refseq     | <a href="#">NM_020780</a>                                                                                                   |
| Protein Refseq  | <a href="#">NP_065831</a>                                                                                                   |
| UniProt ID      | Q9P2K9                                                                                                                      |
| Function        | hedgehog receptor activity; molecular_function                                                                              |