



PTER blocking peptide (CDBP5967)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	PTER phosphotriesterase related [Homo sapiens (human)]
Official Symbol	PTER
Synonyms	PTER; phosphotriesterase related; HPHRP; RPR-1; phosphotriesterase-related protein; parathion hydrolase-related protein; resiniferatoxin-binding, phosphotriesterase-related
Entrez Gene ID	9317
mRNA Refseq	NM_001001484
Protein Refseq	NP_001001484
UniProt ID	Q96BW5
Function	hydrolase activity, acting on ester bonds; molecular_function; zinc ion binding