



RANBP10 blocking peptide (CDBP5978)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RANBP10 RAN binding protein 10 [Homo sapiens (human)]
Official Symbol	RANBP10
Synonyms	RANBP10; RAN binding protein 10; ran-binding protein 10
Entrez Gene ID	57610
mRNA Refseq	NM_020850
Protein Refseq	NP_065901
UniProt ID	Q6VN20
Pathway	Signaling events mediated by Hepatocyte Growth Factor Receptor (c-Met)

Function

Ran GTPase binding; Ran guanyl-nucleotide exchange factor activity; beta-tubulin binding
