



UIMC1 blocking peptide (CDBP5979)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	UIMC1 ubiquitin interaction motif containing 1 [Homo sapiens (human)]
Official Symbol	UIMC1
Synonyms	UIMC1; ubiquitin interaction motif containing 1; RAP80; X2HRIP110; BRCA1-A complex subunit RAP80; receptor associated protein 80; receptor-associated protein 80; retinoid x receptor interacting protein; retinoid X receptor-interacting protein 110; ubiquitin interaction motif-containing protein 1
Entrez Gene ID	51720
mRNA Refseq	NM_001199297
Protein Refseq	NP_001186226
UniProt ID	Q96RL1

Function

K63-linked polyubiquitin binding; histone binding; protein binding; retinoid X receptor binding
