



RASSF10 blocking peptide (CDBP5985)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RASSF10 Ras association (RalGDS/AF-6) domain family (N-terminal) member 10 [Homo sapiens (human)]
Official Symbol	RASSF10
Synonyms	RASSF10; Ras association (RalGDS/AF-6) domain family (N-terminal) member 10; ras association domain-containing protein 10
Entrez Gene ID	644943
mRNA Refseq	NM_001080521
Protein Refseq	NP_001073990
UniProt ID	A6NK89