



ESRP1 blocking peptide (CDBP5988)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ESRP1 is an epithelial cell-type-specific splicing regulator (Warzecha et al., 2009 [PubMed 19285943]).[supplied by OMIM, Aug 2009]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ESRP1 epithelial splicing regulatory protein 1 [Homo sapiens (human)]
Official Symbol	ESRP1
Synonyms	ESRP1; epithelial splicing regulatory protein 1; RBM35A; RMB35A; RNA-binding protein 35A; RNA binding motif protein 35A; RNA-binding motif protein 35A
Entrez Gene ID	54845
mRNA Refseq	NM_001034915
Protein Refseq	NP_001030087

UniProt ID	Q6NXG1
Function	mRNA binding; nucleotide binding