



Rcan2 blocking peptide (CDBP5990)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Rcan2 regulator of calcineurin 2 [Mus musculus (house mouse)]
Official Symbol	Rcan2
Synonyms	Rcan2; regulator of calcineurin 2; Csp2; MCIP2; ZAKI-4; Dscr1l1; calcipressin-2; calcineurin inhibitory protein ZAKI-4; Down syndrome critical region gene 1-like 1; down syndrome candidate region 1-like protein 1; myocyte-enriched calcineurin interactin protein 2; myocyte-enriched calcineurin-interacting protein 2
Entrez Gene ID	53901
mRNA Refseq	NM_001286653
Protein Refseq	NP_001273582
UniProt ID	Q9JHG2

Pathway	Thyroid hormone signaling pathway
Function	nucleotide binding