



REEP3 blocking peptide (CDBP5992)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	REEP3 receptor accessory protein 3 [Homo sapiens (human)]
Official Symbol	REEP3
Synonyms	REEP3; receptor accessory protein 3; C10orf74; receptor expression-enhancing protein 3; receptor expression enhancing protein 3
Entrez Gene ID	221035
mRNA Refseq	NM_001001330
Protein Refseq	NP_001001330
UniProt ID	Q6NUK4