



REEP4 blocking peptide (CDBP5993)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	REEP4 receptor accessory protein 4 [Homo sapiens (human)]
Official Symbol	REEP4
Synonyms	REEP4; receptor accessory protein 4; PP432; C8orf20; receptor expression-enhancing protein 4; receptor expression enhancing protein 4
Entrez Gene ID	80346
mRNA Refseq	NM_025232
Protein Refseq	NP_079508
UniProt ID	Q9H6H4
Function	microtubule binding