



RGP1 blocking peptide (CDBP5995)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RGP1 RGP1 retrograde golgi transport homolog (S. cerevisiae) [Homo sapiens (human)]
Official Symbol	RGP1
Synonyms	RGP1; RGP1 retrograde golgi transport homolog (S. cerevisiae); KIAA0258; retrograde Golgi transport protein RGP1 homolog
Entrez Gene ID	9827
mRNA Refseq	NM_001080496
Protein Refseq	NP_001073965
UniProt ID	Q92546