



RGS22 blocking peptide (CDBP5999)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RGS22 regulator of G-protein signaling 22 [Homo sapiens (human)]
Official Symbol	RGS22
Synonyms	RGS22; regulator of G-protein signaling 22; CT145; PRTD-NY2
Entrez Gene ID	26166
mRNA Refseq	NM_001286692
Protein Refseq	NP_001273621
UniProt ID	Q8NE09
Function	GTPase activator activity