



RHBDD1 blocking peptide (CDBP6000)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RHBDD1 rhomboid domain containing 1 [Homo sapiens (human)]
Official Symbol	RHBDD1
Synonyms	RHBDD1; rhomboid domain containing 1; RRP4; rhomboid-related protein 4; rhomboid-like protein 4; rhomboid domain-containing protein 1
Entrez Gene ID	84236
mRNA Refseq	NM_001167608
Protein Refseq	NP_001161080
UniProt ID	Q8TEB9
Function	endopeptidase activity; protein binding; serine-type endopeptidase activity; serine-type

