



RHBDD2 blocking peptide (CDBP6001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RHBDD2 rhomboid domain containing 2 [Homo sapiens (human)]
Official Symbol	RHBDD2
Synonyms	RHBDD2; rhomboid domain containing 2; NPD007; RHBDL7; rhomboid domain-containing protein 2; H_RG122E10.2a; H_RG122E10.2b; WUGSC:H_RG122E10.2a; WUGSC:H_RG122E10.2b; rhomboid, veinlet-like 7
Entrez Gene ID	57414
mRNA Refseq	NM_001040456
Protein Refseq	NP_001035546
UniProt ID	Q6NTF9

Functionserine-type endopeptidase activity
