



RHBDD3 blocking peptide (CDBP6002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RHBDD3 rhomboid domain containing 3 [Homo sapiens (human)]
Official Symbol	RHBDD3
Synonyms	RHBDD3; rhomboid domain containing 3; PTAG; C22orf3; HS984G1A; rhomboid domain-containing protein 3; pituitary tumor apoptosis
Entrez Gene ID	25807
mRNA Refseq	NM_012265
Protein Refseq	NP_036397
UniProt ID	Q9Y3P4
Function	serine-type endopeptidase activity