



RILPL1 blocking peptide (CDBP6011)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RILPL1 Rab interacting lysosomal protein-like 1 [Homo sapiens (human)]
Official Symbol	RILPL1
Synonyms	RILPL1; Rab interacting lysosomal protein-like 1; RLP1; GOSPEL; RILP-like protein 1; rab-interacting lysosomal-like protein 1
Entrez Gene ID	353116
mRNA Refseq	NM_178314
Protein Refseq	NP_847884
UniProt ID	Q5EBL4