



# RILPL2 blocking peptide (CDBP6012)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

## GENE INFORMATION

|                 |                                                                                                                                                |
|-----------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">RILPL2 Rab interacting lysosomal protein-like 2 [ Homo sapiens (human) ]</a>                                                       |
| Official Symbol | RILPL2                                                                                                                                         |
| Synonyms        | RILPL2; Rab interacting lysosomal protein-like 2; RLP2; RILP-like protein 2; p40phox-binding protein; rab-interacting lysosomal protein-like 2 |
| Entrez Gene ID  | <a href="#">196383</a>                                                                                                                         |
| mRNA Refseq     | <a href="#">NM_145058</a>                                                                                                                      |
| Protein Refseq  | <a href="#">NP_659495</a>                                                                                                                      |
| UniProt ID      | Q969X0                                                                                                                                         |
| Function        | identical protein binding                                                                                                                      |