



RIMS3 blocking peptide (CDBP6015)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RIMS3 regulating synaptic membrane exocytosis 3 [Homo sapiens (human)]
Official Symbol	RIMS3
Synonyms	RIMS3; regulating synaptic membrane exocytosis 3; NIM3; RIM3; regulating synaptic membrane exocytosis protein 3; RIM 3; RIM3 gamma; Rab-3 interacting molecule 3; rab-3-interacting molecule 3
Entrez Gene ID	9783
mRNA Refseq	NM_014747
Protein Refseq	NP_055562
UniProt ID	Q9UJD0

Function	ion channel binding
----------	---------------------
