



# MEX3D blocking peptide (CDBP6018)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |  |
|---------------|--|
| Conjugate     | Unconjugated   |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid   |
| Concentration | 200 µg/mL  |
| Size          | 0.05 mg  |
| Preservative  | None   |
| Storage       | -20°C  |

## GENE INFORMATION

|                 |   |
|-----------------|---|
| Gene Name       | <a href="#">MEX3D mex-3 RNA binding family member D [ Homo sapiens (human) ]</a>  |
| Official Symbol | MEX3D   |
| Synonyms        | MEX3D; mex-3 RNA binding family member D; MEX3; TINO; RKHD1; MEX-3D; RNF193; OK/SW-cl.4; RNA-binding protein MEX3D; RING finger protein 193; bcl-2 ARE RNA binding protein; ring finger and KH domain containing 1; RING finger and KH domain-containing protein 1; ring finger (C3HC4 type) and KH domain containing 1 |
| Entrez Gene ID  | <a href="#">399664</a>  |
| mRNA Refseq     | <a href="#">NM_001174118</a>  |
| Protein Refseq  | <a href="#">NP_001167589</a>  |
| UniProt ID      | Q86XN8  |

**Function**

AU-rich element binding; poly(A) RNA binding; zinc ion binding

---