



MEX3A blocking peptide (CDBP6020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MEX3A mex-3 RNA binding family member A [Homo sapiens (human)]
Official Symbol	MEX3A
Synonyms	MEX3A; mex-3 RNA binding family member A; RKHD4; MEX-3A; RNA-binding protein MEX3A; ring finger and KH domain containing 4; ring finger and KH domain containing protein; RING finger and KH domain-containing protein 4
Entrez Gene ID	92312
mRNA Refseq	NM_001093725
Protein Refseq	NP_001087194
UniProt ID	A1L020

Function	poly(A) RNA binding; zinc ion binding
-----------------	---------------------------------------
