



RPAIN blocking peptide (CDBP6025)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RPAIN RPA interacting protein [Homo sapiens (human)]
Official Symbol	RPAIN
Synonyms	RPAIN; RPA interacting protein; RIP; HRIP; RPA-interacting protein; nuclear transporter; RAP interaction protein
Entrez Gene ID	84268
mRNA Refseq	NM_001033002
Protein Refseq	NP_001028174
UniProt ID	Q86UA6
Function	metal ion binding; protein complex binding