



RSRC1 blocking peptide (CDBP6031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of the serine and arginine rich-related protein family. The encoded protein is involved in both constitutive and alternative mRNA splicing. This gene may be associated with schizophrenia. A pseudogene of this gene is located on chromosome 9. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2012]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RSRC1 arginine/serine-rich coiled-coil 1 [Homo sapiens (human)]
Official Symbol	RSRC1
Synonyms	RSRC1; arginine/serine-rich coiled-coil 1; BM-011; SFRS21; SRp53; serine/Arginine-related protein 53; splicing factor, arginine/serine-rich 21; arginine/serine-rich coiled-coil protein 1
Entrez Gene ID	51319
mRNA Refseq	NM_001271834

Protein Refseq	NP_001258763
UniProt ID	Q96IZ7
Function	protein binding