



## RTN4RL2 blocking peptide (CDBP6032)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

### GENE INFORMATION

|                 |                                                                                                                                                                      |
|-----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">RTN4RL2 reticulon 4 receptor-like 2 [ Homo sapiens (human) ]</a>                                                                                         |
| Official Symbol | RTN4RL2                                                                                                                                                              |
| Synonyms        | RTN4RL2; reticulon 4 receptor-like 2; NgR2; NGRH1; reticulon-4 receptor-like 2; nogo receptor-like 3; Nogo-66 receptor homolog 1; nogo-66 receptor-related protein 2 |
| Entrez Gene ID  | <a href="#">349667</a>                                                                                                                                               |
| mRNA Refseq     | <a href="#">NM_178570</a>                                                                                                                                            |
| Protein Refseq  | <a href="#">NP_848665</a>                                                                                                                                            |
| UniProt ID      | Q86UN3                                                                                                                                                               |
| Function        | receptor activity                                                                                                                                                    |