



RUSC1 blocking peptide (CDBP6036)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RUSC1 RUN and SH3 domain containing 1 [Homo sapiens (human)]
Official Symbol	RUSC1
Synonyms	RUSC1; RUN and SH3 domain containing 1; NESCA; RUN and SH3 domain-containing protein 1; new molecule containing SH3 at the carboxy-terminus
Entrez Gene ID	23623
mRNA Refseq	NM_001105203
Protein Refseq	NP_001098673
UniProt ID	Q9BVN2
Pathway	Trk receptor signaling mediated by the MAPK pathway

Function

SH3/SH2 adaptor activity; actin binding; protein binding
