



RUSC2 blocking peptide (CDBP6038)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a RUN and SH3 domain containing protein that interacts with Rab1b and Rab1-binding protein GM130. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jun 2012]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RUSC2 RUN and SH3 domain containing 2 [Homo sapiens (human)]
Official Symbol	RUSC2
Synonyms	RUSC2; RUN and SH3 domain containing 2; Iporin; iporin; interacting protein of Rab1; RUN and SH3 domain-containing protein 2
Entrez Gene ID	9853
mRNA Refseq	NM_001135999
Protein Refseq	NP_001129471

UniProt ID	Q8N2Y8
Function	Rab GTPase binding; protein binding