



# RWDD3 blocking peptide (CDBP6039)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

## GENE INFORMATION

|                 |                                                                                                             |
|-----------------|-------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">RWDD3 RWD domain containing 3 [ Homo sapiens (human) ]</a>                                      |
| Official Symbol | RWDD3                                                                                                       |
| Synonyms        | RWDD3; RWD domain containing 3; RSUME; RWD domain-containing protein 3; RWD-containing sumoylation enhancer |
| Entrez Gene ID  | <a href="#">25950</a>                                                                                       |
| mRNA Refseq     | <a href="#">NM_001128142</a>                                                                                |
| Protein Refseq  | <a href="#">NP_001121614</a>                                                                                |
| UniProt ID      | Q9Y3V2                                                                                                      |
| Function        | protein binding                                                                                             |