



RWDD3 blocking peptide (CDBP6040)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	RWDD3 RWD domain containing 3 [Homo sapiens (human)]
Official Symbol	RWDD3
Synonyms	RWDD3; RWD domain containing 3; RSUME; RWD domain-containing protein 3; RWD-containing sumoylation enhancer
Entrez Gene ID	25950
mRNA Refseq	NM_001128142
Protein Refseq	NP_001121614
UniProt ID	Q9Y3V2
Function	protein binding