



DLGAP2 blocking peptide (CDBP6045)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The product of this gene is a membrane-associated protein that may play a role in synapse organization and signalling in neuronal cells. This gene is biallelically expressed in the brain, however, only the paternal allele is expressed in the testis (PMID:18055845). Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jun 2014]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DLGAP2 discs, large (Drosophila) homolog-associated protein 2 [Homo sapiens (human)]
Official Symbol	DLGAP2
Synonyms	DLGAP2; discs, large (Drosophila) homolog-associated protein 2; DAP2; SAPAP2; disks large-associated protein 2; PSD-95/SAP90-binding protein 2; discs large-associated protein 2; SAP90/PSD-95-associated protein 2
Entrez Gene ID	9228

mRNA Refseq	NM_001277161
Protein Refseq	NP_001264090
UniProt ID	Q9P1A6
Function	protein binding