



DLGAP3 blocking peptide (CDBP6046)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	DLGAP3 discs. large (Drosophila) homolog-associated protein 3 [Homo sapiens (human)]
Official Symbol	DLGAP3
Synonyms	DLGAP3; discs, large (Drosophila) homolog-associated protein 3; DAP3; SAPAP3; disks large-associated protein 3; DAP-3; PSD-95/SAP90-binding protein 3; SAP90/PSD 95 associated protein 3; SAP90/PSD-95-associated protein 3
Entrez Gene ID	58512
mRNA Refseq	NM_001080418
Protein Refseq	NP_001073887
UniProt ID	O95886

Function	beta-amyloid binding; protein binding
-----------------	---------------------------------------
