



SARM1 blocking peptide (CDBP6049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SARM1 sterile alpha and TIR motif containing 1 [Homo sapiens (human)]
Official Symbol	SARM1
Synonyms	SARM1; sterile alpha and TIR motif containing 1; SARM; SAMD2; MyD88-5; sterile alpha and TIR motif-containing protein 1; tir-1 homolog; SAM domain-containing protein 2; sterile alpha and Armadillo repeat protein; sterile alpha motif domain-containing protein 2; sterile alpha and HEAT/Armadillo motif protein, ortholog of Drosophila
Entrez Gene ID	23098
mRNA Refseq	NM_015077
Protein Refseq	NP_055892
UniProt ID	Q6SZW1

Pathway

Activated TLR4 signalling; Activation of IRF3/IRF7 mediated by TBK1/IKK epsilon; IKK complex recruitment mediated by RIP1; Immune System; Innate Immune System; MyD88-independent cascade; TRAF6 mediated induction of TAK1 complex; TRIF-mediated TLR3/TLR4 signaling
