



# SARS SPIKE blocking peptide (CDBP6058)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

## BACKGROUND

Introduction	The SARS coronavirus, sometimes shortened to SARS-CoV, is the virus that causes severe acute respiratory syndrome (SARS). On April 16, 2004, following the outbreak of SARS in Asia and secondary cases elsewhere in the world.
Keywords	SARS Spike (IN1); SARS; E2; E2 glycoprotein; Human coronavirus spike glycoprotein; Peplomer protein; S; S glycoprotein; Severe acute respiratory syndrome spike glycoprotein; Severe acute respiratory syndrome virus spike glycoprotein; Spike glycoprotein; VGL2; Group IV; Nidovirales; Coronaviridae; Coronavirus; SARS coronavirus; SARS Spike (IN2); SARS Spike (IN3); SARS Spike (CT); SARS Spike Protein; SARS spike glycoprotein; SARS Spike Protein (Middle); SARS (Severe Acute Respiratory Syndrome) Spike protein (Middle);