



## SCARA5 blocking peptide (CDBP6064)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

### GENE INFORMATION

Gene Name	<a href="#">SCARA5 scavenger receptor class A, member 5 [ Homo sapiens (human) ]</a>
Official Symbol	SCARA5
Synonyms	SCARA5; scavenger receptor class A, member 5; Tesr; NET33; scavenger receptor class A member 5; scavenger receptor hlg; testis expressed scavenger receptor; scavenger receptor class A, member 5 (putative)
Entrez Gene ID	<a href="#">286133</a>
mRNA Refseq	<a href="#">NM_173833</a>
Protein Refseq	<a href="#">NP_776194</a>
UniProt ID	Q6ZMJ2

<b>Pathway</b>	Binding and Uptake of Ligands by Scavenger Receptors; Scavenging by Class A Receptors
<b>Function</b>	ferritin receptor activity; scavenger receptor activity