



FBXL20 blocking peptide (CDBP6068)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Members of the F-box protein family, such as FBXL20, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FBXL20 F-box and leucine-rich repeat protein 20 [Homo sapiens (human)]
Official Symbol	FBXL20
Synonyms	FBXL20; F-box and leucine-rich repeat protein 20; Fbl2; Fbl20; F-box/LRR-repeat protein 20; F-box protein FBL2
Entrez Gene ID	84961
mRNA Refseq	NM_001184906

Protein Refseq [NP_001171835](#)

UniProt ID Q96IG2
