



SCRN2 blocking peptide (CDBP6071)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SCRN2 secernin 2 [Homo sapiens (human)]
Official Symbol	SCRN2
Synonyms	SCRN2; secernin 2; Ses2; secernin-2
Entrez Gene ID	90507
mRNA Refseq	NM_001145023
Protein Refseq	NP_001138495
UniProt ID	Q96FV2
Function	dipeptidase activity; molecular_function