



SCRN3 blocking peptide (CDBP6072)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SCRN3 secernin 3 [Homo sapiens (human)]
Official Symbol	SCRN3
Synonyms	SCRN3; secernin 3; SES3; secernin-3
Entrez Gene ID	79634
mRNA Refseq	NM_001193528
Protein Refseq	NP_001180457
UniProt ID	Q0VDG4
Function	dipeptidase activity