



SCUBE3 blocking peptide (CDBP6076)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SCUBE3 signal peptide, CUB domain, EGF-like 3 [Homo sapiens (human)]
Official Symbol	SCUBE3
Synonyms	SCUBE3; signal peptide, CUB domain, EGF-like 3; CEGF3; signal peptide, CUB and EGF-like domain-containing protein 3; CUB and EGF containing protein 3; CUB domain and EGF-like repeat containing 3; signal peptide, CUB and EGF-like domain containing protein 3
Entrez Gene ID	222663
mRNA Refseq	NM_152753
Protein Refseq	NP_689966
UniProt ID	Q8IX30

Function	calcium ion binding; protein binding
-----------------	--------------------------------------
