



SESTD1 blocking peptide (CDBP6088)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SESTD1 SEC14 and spectrin domains 1 [Homo sapiens (human)]
Official Symbol	SESTD1
Synonyms	SESTD1; SEC14 and spectrin domains 1; SOLO; SEC14 domain and spectrin repeat-containing protein 1; protein Solo; huntingtin-interacting protein-like protein
Entrez Gene ID	91404
mRNA Refseq	NM_178123
Protein Refseq	NP_835224
UniProt ID	Q86VW0
Function	NOT 1-phosphatidylinositol binding; phosphatidic acid binding; NOT phosphatidylcholine

binding; phosphatidylinositol-3,4-bisphosphate binding; phosphatidylinositol-3,5-bisphosphate binding; phosphatidylinositol-3-phosphate binding; phosphatidylinositol-4,5-bisphosphate binding; phosphatidylinositol-4-phosphate binding; phosphatidylinositol-5-phosphate binding; NOT phosphatidylserine binding; protein binding
