



SGSM1 blocking peptide (CDBP6089)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SGSM1 small G protein signaling modulator 1 [Homo sapiens (human)]
Official Symbol	SGSM1
Synonyms	SGSM1; small G protein signaling modulator 1; RUTBC2; RUN and TBC1 domain containing 2; RUN and TBC1 domain-containing protein 2; small G protein signaling modulator 1 protein
Entrez Gene ID	129049
mRNA Refseq	NM_001039948
Protein Refseq	NP_001035037
UniProt ID	Q2NKQ1
Function	Rab GTPase activator activity