



SGSM2 blocking peptide (CDBP6090)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SGSM2 small G protein signaling modulator 2 [Homo sapiens (human)]
Official Symbol	SGSM2
Synonyms	SGSM2; small G protein signaling modulator 2; RUTBC1; RUN and TBC1 domain containing 1; RUN and TBC1 domain-containing protein 1; small G protein signaling modulator 2 protein
Entrez Gene ID	9905
mRNA Refseq	NM_001098509
Protein Refseq	NP_001091979
UniProt ID	O43147
Function	Rab GTPase activator activity