



SGSM3 blocking peptide (CDBP6091)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SGSM3 small G protein signaling modulator 3 [Homo sapiens (human)]
Official Symbol	SGSM3
Synonyms	SGSM3; small G protein signaling modulator 3; MAP; RUSC3; RUTBC3; RABGAP5; RabGAP-5; rabGAPLP; merlin binding protein; RUN and SH3 containing 3; merlin-associated protein; RUN and TBC1 domain containing 3; RUN and TBC1 domain-containing protein 3; rab-GTPase-activating protein-like protein; small G protein signaling modulator 3 protein
Entrez Gene ID	27352
mRNA Refseq	NM_015705
Protein Refseq	NP_056520
UniProt ID	Q96HU1

Function

Rab GTPase activator activity; Rab GTPase binding; protein binding
