



SH3RF2 blocking peptide (CDBP6093)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SH3RF2 SH3 domain containing ring finger 2 [Homo sapiens (human)]
Official Symbol	SH3RF2
Synonyms	SH3RF2; SH3 domain containing ring finger 2; HEPP1; POSHER; RNF158; PPP1R39; putative E3 ubiquitin-protein ligase SH3RF2; RING finger protein 158; POSH-eliminating RING protein; SH3 domain-containing RING finger protein 2; heart protein phosphatase 1-binding protein; protein phosphatase 1 regulatory subunit 39; protein phosphatase 1, regulatory subunit 39
Entrez Gene ID	153769
mRNA Refseq	NM_152550
Protein Refseq	NP_689763
UniProt ID	Q8TEC5

Function	ligase activity; phosphatase binding; protein binding; protein phosphatase 1 binding; protein phosphatase inhibitor activity; zinc ion binding
-----------------	--
