



SHISA3 blocking peptide (CDBP6095)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SHISA3 shisa family member 3 [Homo sapiens (human)]
Official Symbol	SHISA3
Synonyms	SHISA3; shisa family member 3; hShisa3; protein shisa-3 homolog; shisa homolog 3
Entrez Gene ID	152573
mRNA Refseq	NM_001080505
Protein Refseq	NP_001073974
UniProt ID	A0PJX4