



SHISA4 blocking peptide (CDBP6096)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SHISA4 shisa family member 4 [Homo sapiens (human)]
Official Symbol	SHISA4
Synonyms	SHISA4; shisa family member 4; TMEM58; C1orf40; protein shisa-4; shisa homolog 4; transmembrane protein 58
Entrez Gene ID	149345
mRNA Refseq	NM_198149
Protein Refseq	NP_937792
UniProt ID	Q96DD7