



SHISA9 blocking peptide (CDBP6097)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SHISA9 shisa family member 9 [Homo sapiens (human)]
Official Symbol	SHISA9
Synonyms	SHISA9; shisa family member 9; CKAMP44; protein shisa-9; shisa homolog 9; cystine-knot AMPAR modulating protein
Entrez Gene ID	729993
mRNA Refseq	NM_001145204
Protein Refseq	NP_001138676
UniProt ID	B4DS77