



SHROOM1 blocking peptide (CDBP6102)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SHROOM family members play diverse roles in the development of the nervous system and other tissues (Hagens et al., 2006 [PubMed 16615870]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SHROOM1 shroom family member 1 [Homo sapiens (human)]
Official Symbol	SHROOM1
Synonyms	SHROOM1; shroom family member 1; APXL2; protein Shroom1; apical protein 2
Entrez Gene ID	134549
mRNA Refseq	NM_001172700
Protein Refseq	NP_001166171
UniProt ID	Q2M3G4

Functionactin filament binding
