



SHROOM3 blocking peptide (CDBP6103)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a PDZ-domain-containing protein that belongs to a family of Shroom-related proteins. This protein may be involved in regulating cell shape in certain tissues. A similar protein in mice is required for proper neurulation. [provided by RefSeq, Jan 2011]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SHROOM3 shroom family member 3 [Homo sapiens (human)]
Official Symbol	SHROOM3
Synonyms	SHROOM3; shroom family member 3; SHRM; APXL3; ShrmL; MSTP013; protein Shroom3; shroom-related protein; F-actin-binding protein
Entrez Gene ID	57619
mRNA Refseq	NM_020859
Protein Refseq	NP_065910

UniProt ID	Q8TF72
Function	actin binding