



SIDT2 blocking peptide (CDBP6104)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SIDT2 SID1 transmembrane family, member 2 [Homo sapiens (human)]
Official Symbol	SIDT2
Synonyms	SIDT2; SID1 transmembrane family, member 2; CGI-40; SID1 transmembrane family member 2
Entrez Gene ID	51092
mRNA Refseq	NM_001040455
Protein Refseq	NP_001035545
UniProt ID	Q8NBJ9
Function	RNA transmembrane transporter activity