



SIGIRR blocking peptide (CDBP6105)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SIGIRR single immunoglobulin and toll-interleukin 1 receptor (TIR) domain [Homo sapiens (human)]
Official Symbol	SIGIRR
Synonyms	SIGIRR; single immunoglobulin and toll-interleukin 1 receptor (TIR) domain; TIR8; single Ig IL-1-related receptor; toll/interleukin-1 receptor 8; single Ig IL-1R-related molecule; single immunoglobulin domain IL1R1 related; single immunoglobulin domain-containing IL1R-related protein
Entrez Gene ID	59307
mRNA Refseq	NM_001135053
Protein Refseq	NP_001128525

UniProt ID	Q6IA17
Pathway	Activated TLR4 signalling; Immune System; Innate Immune System; MyD88:Mal cascade initiated on plasma membrane; Toll Like Receptor 2 (TLR2) Cascade; Toll Like Receptor 4 (TLR4) Cascade; Toll Like Receptor TLR1:TLR2 Cascade; Toll Like Receptor TLR6:TLR2 Cascade
Function	protein binding