



SIGLEC15 blocking peptide (CDBP6107)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SIGLEC15 sialic acid binding Ig-like lectin 15 [Homo sapiens (human)]
Official Symbol	SIGLEC15
Synonyms	SIGLEC15; sialic acid binding Ig-like lectin 15; CD33L3; HsT1361; SIGLEC-15; sialic acid-binding Ig-like lectin 15; CD33 antigen-like 3; CD33 molecule-like 3
Entrez Gene ID	284266
mRNA Refseq	NM_213602
Protein Refseq	NP_998767
UniProt ID	Q6ZMC9
Pathway	DAP12 interactions; Immune System; Innate Immune System